



Anton Paar

PRACTICAL TIPS FOR DENSITY METER
เทคนิคการใช้และการบำรุงรักษาเครื่องวัดความหนาแน่น (DENSITY METER)

11 August 2021

Agenda

- Introduction to Anton Paar (Thailand) Ltd.
- Density principle and instrument
- Practical tips Operating and Cleaning
- Q&A

Anton Paar: History



Margarete Platzer,
Anton Paar's daughter,
launches
**high-precision
manufacturing.**

1932



Presentation of the
first digital density meter
for laboratories (DMA 02C).

1967



Acquisition of: MSB
Morocutti (A), Dr. Wolfgang
Kernchen GmbH (D),
Petrotest GmbH (D), CSM
Instruments (CH).

2007 – 2013



Acquisition
of Quantachrome
Instruments
(Florida, USA).

2018

1922
Locksmith
Anton Paar
establishes the
company in Graz.



1957

Production of the
**first scientific
instrument:** the
Kratky Small-Angle
X-ray Camera.



1989

First steps
in rheometry,
viscometry,
and digestion
technologies.




2016 – 2017

First instruments
for **particle
characterization**,
**Raman
spectroscopy**,
and **atomic force
microscopy.**



Anton Paar: Operating Worldwide



Anton Paar Thailand



Anton Paar (Thailand) Ltd.

**No. 90 CW Tower, 39th Floor,
Unit No. A3902, Ratchadaphisek Road,
HuaiKhwang Sub-District, HuaiKhwang
District, Bangkok 10310**

Office Tel: +66 2 153 9785

Fax: +66 2 153 9786

E-mail: info.th@anton-paar.com

Website: www.anton-paar.co.th

Demo Laboratory





Instruments & Solutions

Business Unit

Measurement



Material characterization



Solution



Business Unit - Measurement

Density & Concentration



Polarimeter & Refractometer



Business Unit - Measurement

Viscometer



Volatile & Consistency



Business Unit - Material characterization

Rheometer



Particle size & Raman



Solid Density & Gas adsorption



Business Unit - Material characterization

Scratch Tester



Atomic Force Microscopes



Coating Thickness



Business Unit - Solution

Process sensor



Microwave Synthesis, Digestion and Extraction



Advantages of the Oscillating U-tube Technique

Measures true density

- No influence of buoyancy in air
- No influence of gravity

Small amounts of sample

- 1 mL
- Easy to control the temperature




Highest precision

- No human influence

Extremely fast

- 1 to 5 min per measurement


For lab and process applications



Density Measurement

- Density determination and adjustment
 - Relation between oscillation period and density
 - Apparatus constants A and B

$$\rho = A \times P^2 - B$$



ρ Density
 P Oscillation period
 A, B Apparatus constants
 (Calculation: using two samples of known density, e.g. air and water)

Power Features

Ahead of Competition



- **U-View™** (crystal clear and high resolution)
- **FillingCheck™** (fast reaction and detection)
- **ThermoBalance™** (best in class fast temperature stability)
- **Full range viscosity correction** (8 x better at highly viscous and dense samples)



launched in
2018

**Pulsed
Excitation
Method**
the ability to perform
measurements on
highly viscous and
dense samples

Checks



Performing a check to ...

... evaluate the accuracy

... detect cleaning problems

Water Check

Air Check

Predefined

Custom Check

Reference Samples

Density (g/cm³)

✓ Passed

0.99810 0.99820 0.99830

0.99819 g/cm³

Density (g/cm³)

✓ Passed


0.00104 0.00114 0.00124

0.00114 g/cm³


Automatically performed with the chosen product settings. A check always refers to a product!

Reports & Results


- Simple pre-defined standard report
- Customize with address line and logo
- GxP compliant if needed




- Check Report
- Adjustment Report
- Measurement Report
- Instrument Settings



Enhanced Data Management with DMA



Practical tips - Standard Operating Procedure



General Info

- Set of written instructions to document a routine within a working process
- Written in a concise, step-by-step, easy-to-read format
 - Improves work consistency, training (personnel changes)

Including

- Was the water/air check passed successfully? If not...
 - Clean and dry the measuring cell
 - Repeat the water check
 - Perform an air/water adjustment if the check failed again
- Has the sample been degassed?
- Was the sample filled without bubbles?
- Is the correct method displayed?


Practical tips - Preventive Maintenance

Before starting a measurement

- Are the adapter connections leaktight?
- Was the water check passed successfully?
- Has the sample been degassed?
- Was the sample filled without bubbles?
- Is the correct method displayed?
- Is the measuring temperature set correctly?

After measurement

- Was the measuring cell cleaned and dried thoroughly?
- Was the air check passed successfully?



The answer to all questions has to be "YES" to get correct results with DMA

Practical tips - Cleaning Recommendations

"Remove the sample from the measuring cell immediately after the measurement, and clean your instrument regularly."



- Make sure that all **wetted parts are resistant** (check in the instruction manual)
- Be aware of **contamination**: Radioactive materials, infectious agents or other materials constituting health hazards
- Samples can chemically react with the previous sample

Filling of Samples

Sample changer

Xsample 320 / 330 / 520 / 530

(Xs 340 / 610 / 630)to come

- Freely programmable sample list
- Plug & Play installation

PFD / SFD

Aerosol Adapter

Heating Attachment



Modularity

Wide range of modules

- Alcolyzer 3001 Beer + Option Color 430nm
- pH 3101/3201, CarboQC ME + Option O₂ Plus, PFD/PFD Plus
- Sample conditioner

- Haze 3001to come



Q&A

www.anton-paar.com